

Basic laboratory techniques

1. Working with glass and automatic pipettes

- a) Weigh an empty glass beaker and write down its mass (m_0) to the table below.
- b) Using the 500 μ l automatic pipette, measure out 5 ml of distilled water into the beaker and weigh it again ($m_{5\text{ ml}}$).
- c) Subtract the weights ($m_{5\text{ ml}} - m_0$) and write the difference ($x_{i=1}$) to the table.
- d) Empty the beaker and repeat the procedure four more times.

Write all the results to the following table:

i	m ₀	m _{5ml}	x _i = m _{5ml} - m ₀	x̄	Δ = (x _i - x̄)	Δ ²
1						
2						
3						
4						
5						

n number of measurements

$$\bar{x} = \frac{\sum x_i}{n} \text{ arithmetic mean}$$

$$s = \sqrt{\frac{\sum \Delta^2}{n-1}} \text{ standard deviation}$$

- e) Repeat the whole procedure using a 1ml glass pipette. Write the results again to the table:

i	m ₀	m _{5ml}	x _i = m _{5ml} - m ₀	x̄	Δ = (x _i - x̄)	Δ ²
1						
2						
3						
4						
5						

Task: Calculate the arithmetic mean and the standard deviation of your weight measurement.